


Routine Maintenance Kit: tube based systems

p/n: RMK-0027



p/n	Description	Quantity
C-QSC10	Quick-fit capillary connectors, pk 10 (connects fused silica to column)	0.2 (2 connectors)
U-DISK1	Filter, disk, sintered PTFE, 5.1 mm section, pk 10	1
U-DISK3	Filter, disk, sintered PTFE, 6.5 mm section, pk 10	1
SERUTD-5065	Split filter tube, stainless steel, 3 1/2 ", packed with charcoal	1
SERUTE-5099	Fused silica transfer line insert (0.25 mm i.d.) and PTFE sleeve, 2m	1
U-FV001	Ferrule, 1/16" × 0.4 mm graphite/Vespel, pk 10 (connects fused silica transfer line)	0.1 (1 ferrule)
U-FV002	Ferrule, 1/16" × 1/16" graphite/Vespel, pk 10 (PEEK connections/bypass line/bulkhead union)	0.6 (6 ferrules)
U-FV003	Ferrule, 1/8 " × 1/16" graphite/Vespel, pk 10 (PEEK connections/bulkhead unions)	0.2 (2 ferrules)
SERZ-0108	Tubing, PEEK, 1/6" x 0.75mm ID, green, 1m	1
U-COV06	O-Ring, low-emission, size 006, pk 10 (Focusing trap)	1
U-COV07	O-Ring, low-emission, size 007, pk 10 (TD autosampler)	1
U-COV10	O-Ring, low-emission, size 010, pk 10 (Tube receivers/caps)	1
U-COV45	O-Ring, low-emission, 5.92 mm i.d., pk 10 (6mm DiffLok Caps)	1